

SIGMA Spectrum Infusion System

Compatible Cleaners

1. 10% Solution of Clorox® Bleach and 90% Water
2. 3M™ HB Quat Disinfectant Cleaner Concentrate
3. A-456 II Disinfectant Cleaner
4. CaviCide ® Disinfectant/Decontaminant Cleaner
5. CaviWipes ® Disinfectant/Decontaminant Cleaner
6. Clorox Commercial Solutions ® Clorox ® Germicidal Wipes
7. Coverage ® Plus Germicidal Surface Wipes
8. DisCide® Ultra Disinfecting Towelettes
9. Dispatch ® Hospital Cleaner Disinfectant Towels with Bleach
10. EnVerros Sanimaster 4®
11. MicroKill Plus Wipes
12. Opti-Cide³ ® Surface Disinfectant Wipes
13. Opti-Cide³ ® Disinfectant and Cleaner Solution
14. Oxivir TB Wipes
15. PDI ® Sani-Cloth® AF Germicidal Disposable Wipe
16. PDI ® Sani-Cloth® Bleach Wipe
17. Sklar Disinfectant Solution Cleaner
18. Super Sani-Cloth ® Germicidal Disposable Wipe
19. Up to 90% Isopropyl Alcohol (IPA)

Note:

- Baxter Healthcare Corporation cannot attest to the safety or effectiveness of the cleaning and disinfection products referenced. The listing is strictly to identify materials compatibility with the SIGMA Spectrum Infusion Pump.
- Consult the Operator's Manual for instruction and cautions for cleaning.

This document is subject to change, please contact Technical Support for additional information:

- Phone: (800) 356-3454 and follow prompts
- Email: MedinaTechSupport@baxter.com

Rx Only. For the safe and proper use of this device, refer to the appropriate operator's manuals.